

NATIONAL HEALTH AND MEDICAL RESEARCH COUNCIL

8.14 RECOMMENDED METHOD OF ANALYSIS FOR THE DETERMINATION OF TOTAL IODINE IN FOODS OTHER THAN MILK

101st session of Council, June 1986

1 SCOPE

This method utilises the catalytic effect of iodide on the destruction of iron (III) thiocyanate by nitrite.

2 APPLICATION

This method is applicable to most foods in general.

3 PRINCIPLE

Ferric ions form a stable orange-red complex with thiocyanate ions. This complex is reduced to a ferrous complex by nitrite ions in the presence of iodide as a catalyst.

The ferrous complex is colourless. At low nanogram quantities of iodide the reduction is proportional to iodide concentration.

4 REAGENTS AND MATERIALS

4.1 General requirements

- 4.1.1 All the reagents should be of analytical grade and distilled water should be used in preference to deionised water.

4.2 Reagents

4.2.1 Standard iodide solutions

- 4.2.1.1 Stock standard iodide 1 g/L:
Dissolve 1.3080 grams potassium iodide in distilled water and dilute to 1 litre in a volumetric flask.
- 4.2.1.2 Working standards:
Dilute appropriately stock standard to obtain working standards of 4, 8, 12, 16 and 20 µg iodide per mL.

- 4.2.2 **Potassium carbonate (30% m/V):**
Dissolve 30 g of potassium carbonate in water and dilute to 100 mL.
- 4.2.3 **Zinc sulphate (10% m/V):**
Dissolve 10 g of zinc sulphate ($\text{ZnSO}_4 \cdot 7\text{H}_2\text{O}$) in water and dilute in 100 mL.
- 4.2.4 **Potassium thiocyanate (0.023% m/V):**
Dissolve 0.23 g of potassium thiocyanate in water and dilute to 1 L.
- 4.2.5 **Sodium nitrite (2.07% m/V):**
Dissolve 2.07 g of sodium nitrite in water and dilute to 100 mL (Stable only for 1 day).
- 4.2.6 **Ammonium iron sulphate reagent:**
Dissolve 77 g of ammonium iron (III) sulphate ($\text{NH}_4 \text{Fe} (\text{SO}_4)_2 \cdot 12\text{H}_2\text{O}$) in approximately 400 mL of water. Add 167 ± 1 mL of nitric acid (S.G 1.42) and dilute to 1 L, warming until all traces of solid dissolve.

5 APPARATUS

- 5.1 Muffle furnace which can be thermostatically controlled to 550°C
- 5.2 Porcelain crucibles, 45mm diameter, with lids.
- 5.3 Centrifuge, with centrifuge tubes
- 5.4 Vortex mixer

6 PROCEDURE

6.1 Dry ashing

- 6.1.1 Into a clean, dry crucible, accurately weigh approximately 1 g of sample with a content of iodine not exceeding 1 ug.
- 6.1.2 Add 1 mL of potassium carbonate solution and then 1 mL of zinc sulphate solution. Slurry the mixture with a glass rod and wash any residue left on the rod back into the crucible.
- 6.1.3 Dry the crucible at 95°C overnight.
- 6.1.4 Place the dried crucible in a muffle furnace at 100°C and cover with a lid.

- 6.1.5 Raise the temperature evenly to 550°C in approximately 90 mins and maintain at this temperature for 1 hour.
- 6.1.6 Remove the crucible and allow to cool to room temperature.
- 6.1.7 Add 1 mL of zinc sulphate solution and slurry the charred residue.
- 6.1.8 Repeat the drying and ashing as before.
- 6.2 Centrifugation
- 6.2.1 Transfer the cooled ash, normally white or grey in colour to a centrifuge tube with 50 ± 0.5 mL of distilled water.
- 6.2.2 Centrifuge for 5 mins to settle insolubles.
- 6.2.3 Decant about half of the clear solution and store in a clean polyethylene container prior to analysis.
- 6.3 Photometric determination
- 6.3.1 Pipette 4 mL of each working standard, containing 4 to 20 ~~mg~~_{ng}/mL iodide, into a 15 x 1.5 cm test tube.
- 6.3.2 A reagent blank is prepared using distilled water.
- 6.3.3 To each tube add 1 mL of distilled water, 1 mL of potassium thiocyanate solution and 2 mL of ammonium iron (III) sulphate solution.
- 6.3.4 Mix the contents of each tube on a vortex mixer.
- 6.3.5 Add, at 90 second intervals, 1 mL sodium nitrite solution and mix again on a vortex mixer.
- 6.3.6 After 20 minutes, at 90 second intervals measure the absorbance of each solution at 450 nm.
- 6.3.7 Sample digestion solutions are treated by the above procedure.
- 6.3.8 It is essential that all solutions should be maintained at the same temperature.
- 6.3.9 All photometric determinations should be carried out in duplicate.

7 CALCULATION AND EXPRESSION OF RESULTS

- 7.1 Prepare a calibration curve from the absorbances of standard solutions and iodide concentration.
- 7.2 From the calibration curve read the iodide concentration of sample solutions and reagent blank.
- 7.3 Calculate the iodide concentration ug/kg, using the following equation.
- 7.4 Iodide Content (ug/kg) = $\left(\frac{(A-B) \times 50}{W} \right)$

Where A = Iodide content of sample solution (^{ug} / mL)

B = Iodide content of reagent blank (^{ug} / mL)

W = Weight of sample (g)

8 REFERENCES

- 8.1 Semi-Automated Method for the Determination of Total Iodine in Food (1980) R.E.D. Moxon and E.J. Dixon: Analyst 105, 344-52.
- 8.2 The Determination of Iodine in Dairy Products (1984) A.M. Chaffey and T.L. Lewis: AGAL Report No 95.